


CHANGE INTERVAL

Recommendation for Changing UROMED Catheter Valves

Catheter valve **with push button**

UROMED Catheter Valve »universal« (REF 1500)
 UROMED Catheter Valve »universal«
 with Catheter Valve Adapter (REF 1501)

change all: **29 days**



UROMED Catheter Valve »compact« (REF 1504)

change all: **29 days**



Catheter valve **with lever**

UROMED Catheter Valve »basic« (REF 1506)

change all: **14 days**



**Please do not hesitate to contact us
 if you have any questions.**

UROMED Kurt Drews KG

Meessen 7/11
 22113 Oststeinbek
 GERMANY

Phone +49 (0) 40 71 30 07-33

Fax +49 (0) 40 71 30 07-99

E-mail homecare@uromed.de

www.uromed.eu



Please find more languages on our website
www.uromed.eu under the heading "Service".